

## Chromatin Immunoprecipitation

Grow cells: 50mL of ~0.5 OD<sub>600</sub> per sample needed (~ 2.5-5 x10<sup>8</sup> cells)

### Day 1:

Collect 50mL of OD<sub>600</sub> 0.5 per sample

For each sample: @RT

Add 1.4 mL formaldehyde (1% final)

Rock for 15min

Add 3 mL 2.5M glycine (150mM final; neutralizes formaldehyde)

Rock for 5 min

(Place on ice if waiting for further time points)

Wash cells:

Spin 5 min @ 3K @ 4°C

Discard supernatant

Resuspend in 25 mL **cold** TBS

Repeat wash

\*\* (the pellet can be snap frozen and kept at -80°C for later)

Resuspend in 400 µl FA-Lysis Buffer (**+ fresh protease inhibitors**)

**(keep samples on ice at all times)**

Prepare eppendorf tubes with 500 µl glass beads each **on ice**

Transfer cells to cold eppendorf tubes

Break cells with bead-beater for 40 min @ 4 °C

Transfer to fresh eppendorf tubes (through a hole in the bottom of the first)

Resuspend any pellet by pipetting

Sonicate:

(Note: all sonicators are different. Times and velocities vary and should be tested for desired result of ~500 bp fragments. These setting are for a Fisher Scientific 550 Sonic Dismembrator)

Sonicate (on #4) for 10 sec ea., transfer to dry ice for 5 sec, then to wet ice ≥ 1 min

Repeat sonication a total of 5x per sample

Add 1mL FA-Lysis Buffer (**+fresh protease inhibitors**)

Spin for 30 min @ max @ 4 °C

Transfer supernatant to new tubes

Spin 1 hr @ max @ 4 °C

Transfer supernatant (the chromatin) to new tubes

\*\* (Snap freeze supernatant in aliquots and store at -80 °C for later)

Split samples:

Inputs: 10 µl + 450 µl TE → store at -20 °C for later

IPs: 100 µl + 300 µl FA-lysis (**+fresh protease inhibitor**) + Ab → O/N rocking @ 4 °C (≥ 4 hrs) [vols. may vary to optimize pulldown]

(Equilibrate Beads overnight, see below)

### **Equilibrate Beads:**

Wash 3x with 1ml TE (spin @10K for 30 sec)

Wash 3x with 1ml FA-lysis buffer

Resuspend to original volume with FA-lysis buffer

### **Day 2:**

IPs only:

Add 6x volumes of Ab ( $\geq 15 \mu\text{l}$ ) of beads (protein-A sepharose if rabbit polyclonal Ab)

Incubate 1.5-2 hrs @ 4 °C on rocker

Wash: @RT

Spin for 30 sec @ 6k xg

Remove supernatant and discard (do **not** disturb bead-bed)

Add 1.5 mL FA-lysis buffer

Rock 5 min [time may vary to optimize]

Repeat washes sequentially with...

2<sup>nd</sup>: FA-500 buffer

3<sup>rd</sup>: LiCl wash buffer

4<sup>th</sup>: TE

Resuspend in 250 $\mu\text{l}$  Elution buffer

Shake for 10 min @ 65 °C

Spin for 2 min @ max

Transfer supernatant to new tube

Add 250 $\mu\text{l}$  TE

IPs and Inputs:

reverse the crosslinking and treat with protease (either order)

Incubate for at least 2 hrs @ 42 °C with 10 $\mu\text{l}$  proteinase K (0.2 ug/ $\mu\text{l}$  final)

Incubate for at least 5 hrs @ 65 °C (usually O/N)

### **Day 3:**

Purify DNA by 2 PCI extractions and 1 CI extraction:

Add equal volume, vortex

Spin 5 min @ max

Extract aqueous layer (top) to a new tube

Transfer 100  $\mu\text{l}$  of input DNA and 400 $\mu\text{l}$  of each sample to new tubes

Add 2.5 volumes EtOH, 1/10 volume 3M NaAcetate, pH5.2, and 1  $\mu\text{l}$  glycogen (20mg/ml)

Mix samples by inversion

Incubate for 30 min @ -80 °C

Spin for 30 min @ max @ 4 °C

Discard supernatant

Dry pellet

Resuspend in water

500  $\mu\text{l}$  for Input DNA

100  $\mu\text{l}$  for IPs

qPCR

## Buffers:

<b>FA-Lysis Buffer</b>	<b>[Final]</b>	<b>[Stock]</b>	<b>1L</b>
Hepes-KOH, pH7.5	50 mM	1 M	50 ml
NaCl	140 mM	5 M	28 ml
EDTA	1 mM	0.5 M	2 ml
Triton X-100	1%	10%	100 ml
Sodium Deoxycholate	0.1%	Powder	1.0 g
Protease Inhibitor cocktail (fresh) (Or Leupeptin & Pepstatin @ 1ug/ml & PMSF @ 1mM)	1X	25X	

<b>FA-500 Buffer</b>	<b>[Final]</b>	<b>[Stock]</b>	<b>500 ml</b>
Hepes-KOH, pH7.5	50 mM	1 M	25 ml
NaCl	500 mM	5 M	50 ml
EDTA	1 mM	0.5 M	1 ml
Triton X-100	1%	10%	50 ml
Sodium Deoxycholate	0.1%	Powder	0.5 g

<b>LiCl Wash Buffer</b>	<b>[Final]</b>	<b>[Stock]</b>	<b>500 ml</b>
Tris-HCl, pH 8.0	10 mM	1 M	5 ml
LiCl	250 mM	5 M	25 ml
NP-40	0.5%	10%	25 ml
Sodium Deoxycholate	0.5%	Powder	2.5 g
EDTA	1 mM	0.5 M	1 ml

<b>2X Elution Buffer</b>	<b>[Final]</b>	<b>[Stock]</b>	<b>50 ml</b>
Tris-HCl, pH 7.5	50 mM	1 M	2.5 ml
EDTA	10 mM	0.5 M	1 ml
SDS	1%	10%	5 ml

Adapted by

Papamichos-Chronakis, M., Petrakis, T., Ktistaki, E., Topalidou, I. & Tzamarias, D. Cti6, a PHD domain protein, bridges the Cyc8-Tup1 corepressor and the SAGA coactivator to overcome repression at GAL1. *Mol. Cell* 9, 1297–1305 (2002). | [Article](#) | [PubMed](#) | [ISI](#) | [ChemPort](#) |

from

M.H. Kuo and C.D. Allis, In vivo cross-linking and immunoprecipitation for studying dynamic Protein: DNA associations in a chromatin environment, *Methods* **19** (1999), pp. 425–433. [Abstract](#) | [PDF \(166 K\)](#) | [View Record in Scopus](#) | [Cited By in Scopus \(277\)](#)